

For products: DM5, DM7

BEFORE USE

- Retain all packaging
- Read the instructions
- Ensure the surface is stable
- Capacity: check ramp label



Scan for more product information

Note: The handle can be moved to the other side of the ramp if this is more convenient. Use a cross-point screwdriver for the fixings and beware of the abrasive upper surface of the ramp. Do not over-tighten.

TO DEPLOY



1. Place the ramp in position with the lip fully on the upper level (see fig 1). The unique design of the end lips holds the ramp in position (see fig. 2).
2. Ensure the ramp does not become dislodged in use. Check that the ramp is firmly located and stable each time it is used.
3. To carry the ramp, use the foldaway carry handle.

MAINTENANCE

If the ramp becomes damaged or slippery it should not be used. The ramp may be cleaned with normal cleaning products. Maximum recommended cleaning temperature: 50°C.

RECYCLING

If you no longer require this product, please liaise with a metals recycling contractor. If this product was loaned to you by a local authority, please contact them to discuss its removal.

WARRANTY

- This product carries a comprehensive 2-year warranty for normal use; conditions apply
- The ramps are not designed for permanent outside use and must not be fixed down
- Use other than published or intended will void the warranty
- Any modification or repair by anyone other than the manufacturer will void the warranty
- Removal of the batch or serial number will void the warranty



Protected by UK patent applications and (European) Registered Community Designs. We undertake a policy of continual review and research, which may lead to changes in product appearance and specification. Date of Issue: April 2025.

Enable Access, Unit R, Cherrycourt Way, Leighton Buzzard, LU7 4UH, UK.

Enable Access Global Ltd, 27 Phibsboro Place, Dublin 7, D07 V20, Republic of Ireland.

For technical support or details on your nearest stockist please call 020 8275 0375.




Code	Max Recom. Height	Length	Capacity	Weight
------	-------------------	--------	----------	--------

Width: 70cm (27.6ins)

DM5	10cm (3.9ins)	56cm (22ins)	350kg (55st)	3.7kg (8.1lbs)
DM7	13cm (5.1ins)	74cm (29ins)	350kg (55st)	4.9kg (10.9lbs)

RAMP LENGTH

a. Measure the step height/level difference.
(Total height of steps or height to resting point on vehicle).







b. Use the **Ramp Length Guide**

Height (up to)	Length		
	1:6	1:8	1:12
10cm (4ins)	60cm (2ft)	80cm (2ft 8ins)	120cm (4ft)
15cm (6ins)	90cm (3ft)	120cm (4ft)	180cm (6ft)
20cm (8ins)	120cm (4ft)	160cm (5ft 4ins)	240cm (8ft)
25cm (10ins)	150cm (5ft)	200cm (6ft 8ins)	300cm (10ft)
30cm (1ft)	180cm (6ft)	240cm (8ft)	360cm (12ft)
45cm (1ft 6ins)	270cm (9ft)	360cm (12ft)	540cm (18ft)
60cm (2ft)	360cm (12ft)	480cm (16ft)	720cm (24ft)

Ramp Length Guide:
Height x Gradient = Length Required

GRADIENT SLOPE

1:6 (10°)*		Acceptable in some cases for short threshold ramps.
1:8 (7°)*		A useful balance between gradient and space suitable for portable ramps in some assisted applications and events, etc.
1:12 (5°)		Often required by social services and local authorities. Minimum required by Building Regulations Document M for installed ramps.
1:15 (4°)		Sometimes required by social services and local authorities for independent access but typically for installed public access ramps.

*Assess safety/risk



Protected by UK patent applications and (European) Registered Community Designs. We undertake a policy of continual review and research, which may lead to changes in product appearance and specification. Date of Issue: April 2025.

Enable Access, Unit R, Cherrycourt Way, Leighton Buzzard, LU7 4UH, UK.

EC REP Enable Access Global Ltd, 27 Phibsboro Place, Dublin 7, D07 V20, Republic of Ireland.

For technical support or details on your nearest stockist please call 020 8275 0375.

