

AKW Advanced Wall
Mounted Extra Wide
Bariatric Fold-up
Moulded Seat w/ Support
Legs Grey Padded Seat
Back & Arms


£793.20 Inc. VAT


Product Code: 04240P


Product Images

















Description

AKW 4000 series extra wide bariatric fold-up shower seat with support legs. The white moulded shower seat has an inset padded seat, backrest and arms, designed to accommodate a variety of individual needs of users, its unique and innovative design provides years of trouble free service. The contrasting grey coloured padding and white mould aids the visually impaired, the padding provides comfort. The seat pad can be easily removed allowing easier cleaning.

At the point of installation, seat height can be adjusted between 390mm – 640mm to maximise comfort by accommodating the clients' needs.

Can be folded away and stored against the wall when not in use.

This AKW shower seat has a maximum load capacity of 254kg / 40 Stone, wall mounted and is designed for easy cleaning.

Key benefits:

- +/- 13mm fine adjustment to support leg, allowing for gradient in floor to optimise comfort on a level seating
- Seat pad easily removable for cleaning
- Folds compactly against the wall when not in use
- Ideal for situations where side transfers are required
- Independent folding arms makes it suitable for side transfer from either side
- Ribbed design on arms and seat pads prevents cushion from moving during side transfer
- Cushion fixing lugs prevent the cushion from moving during side transfer

Technical details:

- Maximum load capacity (kg / Stone): 254 / 40
- Armrest colour: Grey
- Seat colour: Grey
- Back colour: Grey
- Projection when folded (mm): 250
- Projection of seat from wall (unfolded) (mm): 575
- Width of seat (mm): 580
- Maximum height (mm): 640
- Minimum height (mm): 390
- Materials: Rust-proof stainless steel

Warranty:

• Lifetime manufacturer warranty

Note:

- The lifetime warranty excludes normal wear and tear in general domestic applications and in accordance with AKW installation instructions. It applies only to products being installed for the first time and according to relevant local and national regulations
- It is the responsibility of the installer to install the seat in accordance with best practice
- Support legs must not exceed a maximum of 1 in 40 shower base slopes when placed in the down position. Curved surfaces such as corner of bathtubs are not acceptable
- Wall fixings are not included

Additional Information

| Manufacturer | AKW / Contour |
|------------------------|-----------------------|
| Operation | Manual Folding |
| Configuration | Seat, Arms & Backrest |
| Seat Width (mm) | 500 - 599 |
| Height Adjustment (mm) | 0 (Fixed) |
| Folding Seat | Yes |
| Padded | Yes |