






Independent 4 Life


AKW Advanced Wall Mounted Extra Wide Fold-up Moulded Horseshoe Blue Padded Seat w/ Support Legs

£318.71 Inc. VAT

Product Code: 04120P

Product Images





Description

AKW 4000 series standard extra wide horseshoe fold-up shower seat with support legs. The white moulded shower seat has an inset padded horseshoe seat, extra wide design accommodates for a variety of individual needs making it easier for intimate hygiene. The contrasting blue coloured padding aids the visually impaired. The seat pad can be easily removed making cleaning easier.

At the point of installation, seat height can be adjusted between 390mm – 640mm to maximise comfort by accommodating the clients' needs.

Can be folded away and stored against the wall when not in use, allowing the shower to be accessible to all.

This AKW shower seat has a maximum load capacity of 254kg / 40 Stone, wall mounted and is designed for easy cleaning.

Key benefits:

- Made from stainless steel, ensuring durability for long-term use
- Comfortable in use
- +/- 13mm fine adjustment to support leg, allowing for gradient in floor to optimise comfort on a level seating
- Seat pad is easily removeable for easy cleaning
- Folds compactly against the wall when not in use
- Ideal in situations where side transfer is required

Technical details:

- Maximum load capacity (kg / Stone): 254 / 40
- Seat pad colour: Blue
- Projection when folded (mm): 155
- Projection of seat from wall (unfolded) (mm): 486
- Width of seat (mm): 578
- Maximum height (mm): 640

www.independent4life.co.uk | info@independent4life.co.uk | Tel: 0800 500 3001 | VAT No: 942447711

EORI GB942447711000

- Minimum height (mm): 390
- Materials: Rust-proof stainless steel

Warranty:

- Lifetime manufacturer warranty

Note:

- The lifetime warranty excludes normal wear and tear in general domestic applications and in accordance with AKW installation instructions. It applies only to products being installed for the first time and according to relevant local and national regulations
- It is the responsibility of the installer to install the seat in accordance with best practice
- Support legs must not exceed a maximum of 1 in 40 shower base slope when placed in the down position. Curved surfaces such as corner of bathtubs are not acceptable
- Wall fixings are not included

Additional Information

Manufacturer	AKW / Contour
Operation	Manual Folding
Configuration	Seat Only
Seat Width (mm)	500 - 599
Height Adjustment (mm)	0 (Fixed)
Folding Seat	Yes
Padded	Yes