

AKW Advanced Wall


Mounted Extra Wide Fold-


up Moulded Horseshoe, w/ Support Legs - Grey Padded Seat, Arms & Back

Product Code: 04250P

Product Images







Description

AKW 4000 series standard horseshoe fold-up shower seat with support legs. The white moulded shower seat has an inset grey padded horseshoe seat, backrest and arms, designed for accommodating varied individual needs, making it easier for intimate hygiene. The contrasting Grey padding and white mould provides aid for visual impairment, the padding provides comfort. the seat pad can be easily removed making cleaning easier.

At the point of installation, seat height can be adjusted between 390mm – 640mm to maximise comfort by accommodating the clients' needs.

Can be folded away and stored against the wall when not in use, allowing the shower to be accessible to all.

This AKW shower seat has a maximum load capacity of 254kg / 40 Stone, wall mounted and is designed for easy cleaning.

Key benefits:

- Made from stainless steel, ensuring durability for long-term use
- Comfortable in use
- Extra wide, accommodating a variety of individuals
- Horseshoe shape accommodates a variety of needs
- +/- 13mm fine adjustment to support leg, allowing for gradient in floor to optimise comfort on a level seating
- Seat pad is easily removable for cleaning
- Folds compactly against the wall when not in use
- Independent folding arms makes it suitable for side transfer from either side
- Ribbed design on arms and seat pads prevents cushion from moving during side transfer

Technical details:

- Maximum load capacity (kg / Stone): 254 / 40
- Armrest colour: Grey
- Seat colour: Grey

- Back colour: Grey
- Projection when folded (mm): 255
- Projection of seat from wall (unfolded) (mm): 575
- Width of seat (mm): 580Maximum height (mm): 640
- Minimum height (mm): 390
- Materials: Rust-proof stainless steel

Warranty:

Lifetime manufacturer warranty

Note:

- The lifetime warranty excludes normal wear and tear in general domestic applications and in accordance with AKW installation instructions. It applies only to products being installed for the first time and according to relevant local and national regulations
- It is the responsibility of the installer to install the seat in accordance with best practice
- Support legs must not exceed a maximum of 1 in 40 shower base slope when placed in the down position. Curved surfaces such as corner of bathtubs are not acceptable
- Wall fixings are not included

Additional Information

Manufacturer	AKW / Contour
Operation	Manual Folding
Configuration	Seat, Arms & Backrest
Seat Width (mm)	500 - 599
Height Adjustment (mm)	0 (Fixed)
Folding Seat	Yes
Padded	Yes